

# PRODUCTION OF PARTICULATE AND ELECTROMAGNETIC RADIATION

---

Shiva Das, Ph.D.  
Duke University

# BASIC CONCEPTS OF X-RAYS

---

# Three Types of Photon Productions

- Radioactive Decay:  $\gamma$  ray with a fixed energy.
- Characteristic X-ray: knock out an orbital electron to produce characteristic x-ray with a fixed energy.
- Bremsstrahlung

# Bremsstrahlung X-ray

- When a high speed electron passes near the nucleus, it suffers a sudden deflection and acceleration due to the attractive Coulomb force between the electron and the nucleus.
- A part or all of the electron energy is dissociated from the electron and propagates as radiation.
- Multiple bremsstrahlung interactions may happen to an electron.
- The resulting bremsstrahlung photon may have any energy up to initial electron energy (not a fixed energy!).


# Bremsstrahlung X-ray

- The direction of bremsstrahlung photons depends on the energy of incident electrons
- When  $e^-$  energy  $\sim 100$  keV, bremsstrahlung x-rays are emitted almost in all directions
- As the  $e^-$  energy increases, the direction of x-rays becomes increasingly forward
- At MeV  $e^-$  energy, most photons are in the forward direction
- A transmission target is used in Linear Accelerator

# X-ray Energy Spectra

- X-rays are heterogenous in energy
  - A continuous distribution of energies for X-rays + discrete characteristic x-rays
- Filtration effect
  - inherent filtration: remove the low energy x-rays
  - added filtration: enhance the penetration power
- Average X-ray energy
  - $\frac{1}{3}$  of the maximum energy

## Low Energy Therapy Units



### “Grenz-Ray” Therapy

Described as soft or low energy x-rays


X-ray potentials of 20 kV or less

Comparable to mammography energies

No longer used due to extremely low penetration

HVL  $\approx$  0.04 mm Al

## Low Energy Therapy Units



### Contact therapy

Operates at potentials of 40-50 kV

Soft component of beam absorbed by .5-1.0 mm thick aluminum filtration

Used with SSD of approximately 2cm

Provides a rapid decrease in depth dose

Incident beam has maximum dose at skin surface and underlying skin is “spared”

HVL  $\approx$  1.5 mm Al





Fig. 15-7 A handheld contact-therapy machine used to treat superficial skin lesions. The operators, one to monitor the patient and the other to hold the applicator, must wear protective shielding during the treatment application.



## Low Energy Therapy Units



### Superficial therapy

Operating potential of 50-150 kV

Filtration of 1-6 mm aluminum to harden beam


Treatments usually given with cone applicators  
of 15-20 SSD

Limited to treating tumors of  $\sim 5$  mm depth  
(90% depth dose)

HVL  $\approx 1-8$  mm Al



# Low Energy Therapy Units



## Orthovoltage therapy


Operating potentials of 150-500 kV

Most equipment operated in ranges of  
200-300 kV and 10-20 mA

SSD usually set at 50 cm and has  
adjustable field size

Maximum dose close to surface with 90%  
dose line at 2 cm depth


HVL  $\approx$  1-4 mm Cu



## Low Energy Therapy Units




### Depth Dose Characteristics



- a) Grenz rays: 0.04 mm Al
- b) Contact therapy: 1.5 mm Al
- c) Superficial: 1-8 mm Al
- d) Orthovoltage: 1-4 mm Cu
- e)  $^{60}\text{Co}$ : 1.2 cm Pb

## Megavoltage Therapy Units



$\gamma$ - ray emitters from radionuclides as they undergo radioactive disintegration, include:

$^{226}\text{Ra}$

$^{137}\text{Cs}$

$^{60}\text{Co}$

$^{60}\text{Co}$  most suitable because:

Higher possible specific activity (curies per gram)

Greater radiation output per curie

Higher average photon energy

## Cobalt Teletherapy Units


Produced in nuclear reactor by irradiating stable  $^{59}\text{Co}$  with neutrons

$^{60}\text{Co}$  usually in forms of cylinders, discs or pellets


A double welded seal prevents  $\beta$  particle leakage (.32 MeV)

Decays in 5.26 y to stable  $^{60}\text{Ni}$

Gamma rays of 1.17 & 1.33 MeV used in therapy



## Cobalt Teletherapy Units



El Dorado Unit



Theratron Unit



## Cobalt Teletherapy Units



Theratron  
780



Theratron  
1000



Radon Medical,  
Turkey

# Accelerators

- Accelerators are initially built to accelerate charged particles for nuclear and high energy physics research.
- There are two types of accelerators:
  - Straight line accelerators
  - Cyclic accelerators



## Straight line

## Cyclic

<b>Linac: electrons</b>	<b>Synchrotron: electrons, protons</b>
<b>X-ray tube: electrons</b>	<b>Betatron: electrons</b>
	<b>Cyclotron: Proton, ions, deuterons</b>
	<b>Microtron: electrons</b>

# Linear Accelerators


- Medical linear accelerators accelerate electrons with a kinetic energies from 4 to 25 MeV by microwave power source.
- The microwave power source is usually defined as electromagnetic energy ranging from approximately 1 GHz to 100 GHz in frequency.

# Linear Accelerators

## Five basic Components


- Injection system – electron source
- RF system- magnetron/klystron, modulator, modulator control unit, waveguide & circulator
- Auxiliary system – vacuum pumping system (waveguide), water cooling system, gas dielectric system for transmitting microwaves from RF generator to accelerating waveguide.
- Beam transport system – target/scattering foil, magnetic steering, focusing
- Beam collimation & monitoring system – shaping, monitoring of beams

# Linear Accelerator




A block diagram of Linear accelerator

# Simple LINAC




# Modern LINAC



# Microwave Generators

- Magnetron and Klystron are two typical microwave generators normally used in Linear accelerators
- Magnetron:
  - Supply 2 MW peak power, more compact, and inexpensive
  - Can be mounted in the Gantry
- Klystron:
  - Supply higher peak power (5 MW), bulky, and more expensive
  - require a low RF power as input, Usually installed in or behind the gantry stand


# Magnetron



- Electrons spiral outward from cathode in vacuum.
- Induce high power RF fields in cavities.
- Converts up to 60% to microwave power, extracted with an output antenna.
- High heat limits microwave power.



# Klystron



- Electrons are accelerated/decelerated by oscillating RF field across buncher cavity, causing bunching of electrons.
- Catcher cavity (same resonant frequency as buncher cavity) is excited by electrons and create high power RF field in catcher cavity.
- Bremsstrahlung created as a result of collision with collector.

## **Magnetron vs. Klystron:**


- Magnetron is smaller
- Magnetron does not require input RF signals
- Klystron mounted within a tank of insulating oil (has to be stationary, cannot be within rotating gantry).
- Magnetron less stable.

## **Klystron:**

- Klystron needs to produce MW power – not possible to continuously operate.
- Operates at a duty cycle (fractional on-time) of  $10^{-3}$  to  $10^{-4}$ .
- Average power KW, instantaneous power MW.


# What Is a Waveguide?

- A waveguide is gas filled (if only used for EM waves) or evacuated (if EM + electron acceleration) metallic structure of rectangular or circular cross section
- Gas filled waveguide from Klystron/Magnetron transmits microwave RF source into accelerating waveguide.
- Inside the evacuated accelerating waveguide, electrons are accelerated to a near light speed.



## Pulsed modulator

- Housed in the linear accelerator vault.
- Pulse lengths of a few  $\mu\text{s}$  ( $\sim 1000$  pulses per second)
- Pulse transformer is used to amplify voltage
- Pulses sent to the cathode of Klystron/Magnetron as well as cathode of electron gun.




electron gun


# Wave and Particle Duality

- Waveguide theory is an example of how particle (electrons) and electromagnetic field (Microwave) exchange energy.
- The propagation of microwaves through a uniform waveguide is governed by Maxwell's equations and boundary conditions.

(a)




(b)



# Traveling and Standing Wave-guides

- Traveling wave-guide:
  - The residual power is absorbed by a dummy load
  - No reflected wave
- Standing wave-guide:
  - A special device (isolator/circulator) to deal with unwanted reflected power going back into Klystron.
  - Shorter wave-guide can be used to achieve the same electron energy as the traveling waveguide.
  - Requires 25% greater RF power








# Beam Delivery


- Bending Magnet
- Target
- Fixed (primary collimator)
- Flattening filter
- Monitor ionization chambers
  - Transmission type
- Movable collimators and Multi-leaf collimators



## Bending magnet

- For energies  $> 6$  MV, waveguide is too long and must be placed horizontal, so beam has to be bent.
- $90^\circ$ ,  $270^\circ$ ,  $112.5^\circ$  bending magnet.
  - $90^\circ$  magnet – acts as a spectrometer that bends higher energies less, so large focal spot;
  - $270^\circ$  magnet – refocuses spectral spread to provide a small focal spot, but bulky (height);
  - $112.5^\circ$  magnet – advantages of small space with small focal spot.

# Electron Transport System



# Target

Electrons strike target to create: (1) collisional losses; (2) radiative losses; (3) scattering.


- Collisional losses result in heat – dissipated using circulating chilled water.
- Radiative losses produce photons that are in the same direction as the incident electrons.
- Scattering of electrons in the target result in angular spread of the resulting radiation. The higher the target atomic number, the higher the scattering.

## Target selection:

- Intermediate to low Z targets produce higher photon yield, but need to be thicker to absorb all electrons (e.g., lead 10mm, tungsten 5 mm, aluminum 5 cm – too thick to accommodate in head).
- Copper is intermediate Z and suitable.



# Flattening filter

- Flattening filter converts the gaussian shaped field profile into a “flat” beam profile
- 40 x 40 field (100 SSD) at 10 cm depth should have dose variation within the central 80% of field not more than  $\pm 3\%$  of central axis value.
- Increased hardening in central portion results in “horns” at depths less than 10 cm (horns at  $d_{max} \pm 6\%$ )



# Clinical electron beams

- Target and flattening filter are moved out.
  - Beam output is a pencil beam with gaussian profile – cross sectional area too small for clinical use.
  - Scattering foil on rotating carrousel is used to broaden beam. Thin foil of high Z material, e.g., copper or lead.
    - Thin to prevent secondary photons.
    - Resulting beam is still gaussian, but broader.
  - Electron applicator, in conjunction with jaw position is used to further flatten beam.
- 
- ❖ Secondary photons are caused by interaction of electrons with scattering foil, ion chamber, jaws and electron applicator.
  - ❖ Scattering foil, ion chamber and applicator should be as low Z as possible to reduce X-ray contamination. X-ray contamination can range from 1% for 6 MeV to 5% for 25 MeV.







# Multi-leaf Collimator Designs

- Each manufacturer has a different design of MLC
  - Location, leaf width, and leaf end design
  - Single focused or double focused
  - Restrictions on motion (path, over-travel, interleaf)
  - Field size
- These factors have an impact on IMRT delivery and must be considered in treatment planning

# Double Focused MLC




**Focused in in-plane (Y)**




**Focused in cross plane (X)**

# Single Focused MLC



**Focused in in-plane ( Y )**



**Focused in cross-plane ( X )**

# Beam Collimation and Penumbra


- Transmission Penumbra
  - The inner surface of the collimator is parallel to the CAX
- Geometric Penumbra
  - Due to the finite radiation source , s

$$P_d = s (SSD + d - SDD) / SDD$$


Where SDD is the source to diaphragm distance

- Physical Penumbra
  - define the distance between the 20%-90% IDL

# Rounded Leaf End vs Penumbra




# Siemens MLC




# Varian MLC System

Varian  
Tertiary  
MLC



# Elekta MLC System






# Physical Leaf Length vs. Over-travel Distance


- The MLC physics leaf length (project to iso-center) is 16 cm, 30 cm, 32.5 cm for Varian, Siemens, and Elekta Accelerators, respectively.
- The distance that each individual leaf passes over iso-center is called over-travel distance, without leaving an uncovered region behind the leaf.

# Over-travel Distances

- For Siemens and Elekta machines, the over-travel distances are 10 cm and 12.5 cm, respectively, without leaving a uncovered region behind the leaf.
- For Varian MLC, X jaw is used to cover the uncovered region of MLCs. The over-travel distance = 2 cm (x jaw over-travel distance) + 15 cm = 17 cm.
- The maximum differences between the leading leaf and the trailing leaf  $\leq 15$  cm

# Tongue & Groove






a




b




c


# Leaf Motion Constraints



**Interdigitation (Varian)**



**No interdigitation  
(Siemens)**



**Minimum Gap (Elekta)**

**Segmentation is  
affected by these  
constraints**

# MLC Leakages




- Intra-leaf leakage
  - Thickness (1-2%)



- Inter-leaf leakage
  - Tongue & groove (2-4%)
- Leaf-end leakage
  - Rounded leaf end (20%)
  - Flat leaf end (1-2%)





**Intra-leaf leakage  
0.8%**

**Inter-leaf leakage  
(1.5%)**


**Leaf end leakage  
(1.5%)**

**Siemens Primus  
6MV**

# MLC Leakage and Backup Jaws

- MLC leakage can be minimized by letting backup jaws following each IMRT segment.
  - Varian: Backup jaws do not follow each MLC segment.
  - Siemens: Backup jaws follow each segment.
  - Elekta: Backup Jaws follow each segment.







Huq, MS, et.al. *Phys Med Biol* 2002; 47: N159-70.

# PROTON THERAPY

- Proton beams exhibit a Bragg peak effect, depositing a relatively large proportion of their initial energy at the end of their range.
  - Bragg peak, in practice, is blurred by scattering foils, multiple coulombic scattering.
  - Small fields exhibit a lower ratio of Bragg peak to entrance dose due to lateral side scatter.
- 
- Proton beams have very sharp lateral penumbra for low energy protons and high energy protons (shallow depth).
  - Allows for sharp separation between beam path and surrounding OARS - Big factor in deciding to treat patients with protons (aside from Bragg peak).



# Requirements

- Energy of 200 – 250 MeV for penetration (e.g., prostate bilateral prostate treatment can be achieved with energies of 225 MeV).
- Supply about  $5 \times 10^{11}$  protons/minute for an average beam delivery time of less than 3 min per patient.


# Equipment

- Duoplasmatron source
- Radio frequency quadrupole (RFQ) accelerator
- Synchrotron
- Beam line to treatment rooms

## Duoplasmatron source


- Hydrogen gas is injected into high pressure chamber
- Filament creates electrons that mix with hydrogen to create a voltage arc that heats the hydrogen and turns it into a plasma.
- Solenoidal magnet is used to focus the electrons and protons created by plasma to pass through small opening into second chamber.
- Plasma cools in second chamber and electrons are attracted to chamber walls.
- Protons exit plasma surface and are accelerated by high voltage electrodes to RFQ accelerator.

# Duoplasmatron source




# Radiofrequency quadrupole (RFQ) accelerator

- Enclosed box with 4 pole pieces having ridges down entire length.
- Radiofrequency power applied to poles produce axial electric field that accelerates protons and transverse magnetic field for focus.
- Proton energy at exit is about 2 MeV.







Protons exiting RFQ accelerator are turned  $180^\circ$  with bending magnet and injected into synchrotron.


# Synchrotron

- Eight dipole magnets (2 at 4 corners) bend the protons around the synchrotron.
- RF cavity accelerates the protons at each pass.
- Protons gain energy with each pass because of RF acceleration, so with each pass the dipole magnet strength must increase and RF frequency must increase.
- To reach energy of 250 MeV requires about 2.8 million passes.
- Once desired energy is reached, dipole magnetic strength and RF frequency are held constant to maintain energy.
- Beam is extracted to switchyard and then to treatment rooms.



# Proton therapy switchyard

Treatment rooms could have rotating gantry or fixed gantry (fixed vertical, fixed horizontal).



# Methods to produce laterally uniform fields

- Proton beam that comes out of synchrotron is a few mm wide.
- Need to be widened and made uniform across field for clinical use.


## Techniques

- Scattering foils
- Raster scanning

# Methods to produce laterally uniform fields


## **Scattering foils:**

- Lead scatterer can be used to widen beam, but results in a gaussian (not flat) profile.
- A second scatterer that is approximately gaussian in shape is used to flatten beam.



vacuum pipe


pre-collimator




1st scatterer  
wedges



gaussian  
fluence profile



2nd scatterer  
carousel




flattened  
fluence profile

# Methods to produce laterally uniform fields

## Raster scanning:

Narrow beam is steered across the field using magnets.






# Modulating penetration of beam

- Tumors are at different depths in patients, so beam energy must be adjusted on a patient-specific basis.
- Within a specific patient, tumor maximum depth can vary across the field, so beam energy must be adjusted laterally.
- Within a specific patient, tumor minimum-to-maximum depth can vary across the field, so beam energy must be dynamically varied even at every point in the field.

# Modulating penetration of beam

*Tumors are at different depths in patients, so beam energy must be adjusted on a patient-specific basis.*

- Range shifter spiral wedge can be used to tune the energy for a specific patient.
- Adds scatter and increases energy spread in the proton beam (creates undesirable energies).
- A bending magnet with collimator slits can be used to remove the undesirable energies.



# Modulating penetration of beam


*Within a specific patient, tumor maximum depth can vary across the field, so beam energy must be adjusted laterally.*

- A plastic bolus with varying depths across the field using a computer-controlled milling machine that is linked to the treatment planning system.

# Modulating penetration of beam

*Within a specific patient, tumor minimum-to-maximum depth can vary across the field, so beam energy must be dynamically varied even at every point in the field.*

- Spinning “propeller” inserts different thicknesses in front of beam, smearing the energy.




## Final collimation at patient

- Similar to electrons – collimation needs to be done as close to patient as possible because of multiple coulombic scattering in bolus and air.
- Multiple cones are available – custom apertures are placed in cone to define field shape.
  - Cast from Lipowitz metal poured in a hotwire cut styrofoam mold.
  - For high precision, milled precisely from soft metal such as brass.

# Intensity Modulated Proton Therapy


- Can be achieved using a range shifter, bolus (not necessarily needed), propeller and some method to control the beam fluence across the field.
- Beam fluence control.
  - Using an MLC.
  - Magnetic scanning of the beam in a controlled fashion to vary the fluence across field.



# Other accelerators

## Cyclotrons

- Acceleration of ions (protons, deuterons, heavier ions)
- Two evacuated Dees.
- Uniform magnetic field.
- Alternating electric field between electrodes that accelerate particles every time they cross electrodes.
- Magnetic field causes them to bend around.
- As particle energy increases, radius of bending increases, causing a spiral pattern to be traced.
- Can accelerate protons to about 10 MeV and deuterons to about 20 MeV.



Precision Graphics

## Other Accelerators: Cyclotrons



Cyclotrons are used for: isotope production, particle therapy, neutron therapy (by hitting protons or deuterons on a beryllium target)

